

Preparing standard solutions

1a Direct weighing method


REQUIREMENTS FOR PRIMARY STANDARD

- Soluble in water
- X react w/ air
- X hygroscopic
- non-volatile
- non-toxic

cannot be performed for	$I_2(s)$	volatile \rightarrow sublimates
	$KMnO_4(s)$	Strong OA, reactive
	$NaOH(s)$	neutralises by reacting w/ CO_2 hygroscopic


STEPS

- Eg. prepare 0.1M NaCl solution



weigh 1.465g of solid accurately w/ electronic balance


molarity \cdot volume = mole $0.1 \cdot 0.25 = 0.025 \text{ mol}$ \rightarrow volumetric flask capacity
mole \cdot Mr = mass $0.025 \times (23.1 + 35.5) = 1.465 \text{ g}$
 \swarrow Na \searrow Cl



Dissolve NaCl(s) in water completely.
(如果溶不到用 stirrer / 加热)

为啥不直接倒进 volumetric flask 里?

- flask 的开口太小, NaCl 会粘在开口上
- 倒进 flask 以后很难让不是 very soluble 的物质落在水里




Pour solution in beaker to the volumetric flask.

Rinse the beaker w/ distilled water.

Transfer all the washings to the volumetric flask.

避免还有 NaCl 遗在了 beaker 里



Add distilled water to volumetric flask until it reaches 250 cm^3 graduation mark. Add a stopper then Shake well to dissolve.


1b Double standardisation

OBJECTIVES & METHOD

- find actual concentration of std. solⁿ
- eg. After weighing, some Na_2CO_3 reacts w/ CO_2 .
 \hookrightarrow molarity of resultant solⁿ will be smaller than original.
- find by titrating against other std. solⁿ.

STEPS

- eg. prepared 1.0M Na_2CO_3 solution \rightarrow titrate against standard 1.0M HCl



used avg. of 48.6 cm^3 1.0M HCl to change colour (red \rightarrow orange)

$HCl : Na_2CO_3 = 2 : 1$
$1.0 \times 0.0495 : x \times 0.025 = 2 : 1$
$x = 0.972$

2 Diluting method

STEPS

- dilute solⁿ of known molarity to specific volume.
- 限制: 只可以使溶液变稀, 不可以变浓
- eg. prepare 100 cm^3 of 0.01M $Na_2CO_3(aq)$ w/ 250 cm^3 of 0.1M $Na_2CO_3(aq)$.
 1. Using a 10.0 cm^3 pipette, transfer 10 cm^3 of 0.1M $Na_2CO_3(aq)$ into 100.0 cm^3 volumetric flask.
 2. Add distilled water into the volumetric flask until it reaches graduation mark.
 3. Stopper the flask and invert it several times to mix the contents well.

PRACTICE

- 0.1M 250.0 cm^3 solⁿ \rightarrow 0.01M 250.0 cm^3 solⁿ
 1. Using 25.0 cm^3 pipette, transfer 25 cm^3 of 0.1M solⁿ into 250 cm^3 volumetric flask.
 2. Add distilled water into volumetric flask until graduation mark is reached.
 3. Stopper the flask, invert it several times to mix the content well.
- 0.1M 250.0 cm^3 solⁿ \rightarrow 0.02M 100.0 cm^3 solⁿ
 1. Using 20.0 cm^3 pipette, transfer 20 cm^3 of 0.1M solⁿ into 100 cm^3 volumetric flask. \rightarrow 为什么不用两次 10.0 cm^3 pipette?
 - random error \uparrow
 - \hookrightarrow 人为误差 (vs systematic error 系统误差)
 - \hookrightarrow burette reading having $\pm 0.05 \text{ cm}^3$ max. error
 2. Add distilled water into volumetric flask until graduation mark is reached.
 3. Stopper the flask, invert it several times to mix the content well.