
# Titration (calculations)

## 1 Questions

#### LEVEL I

25.0 cm³ CH<sub>3</sub>COOH is withdrawn to 250.0 cm³ v-flask. Sol. in v-flask is withdrawn by 25.0 cm³ pipette to conical flask, and titrated against O.IM NaOH. 30 cm³ NaOH is needed for complete reaction. Find Stock CH<sub>3</sub>COOH concentration.



$$CH_{3}COOH n: NOOH n = 1:1 = n: 0.003$$


$$N = 0.003 \xrightarrow{CH_{3}COOH used in}$$

$$\frac{x}{10} = \frac{0.003}{0.025}$$

$$x = 1.2 \quad mol \quad dm^{-3}$$

## LEVEL 1

0.6 g monobasic acid is dissolved completely in a beaker W/100 cm³ distilled water. The sol. in beaker is then poured into 250.0 cm³ v-flask. Distilled water is added until the graduation mark is reached. 10.0 cm³ of solution in v-flask is pipetted into conical flask. 4 cm³ o.1 M NaoH is required for complete reaction. Find acid Mr.



母数
$$ACID N: NOOH N = 1:1 = \frac{n}{25}:0.0004$$

$$N = 0.01$$

$$\frac{0.69}{Mr} = 0.01$$

$$Mr = 60.09 mol^{-1}$$

### LEVEL III-

0.75 g of Solid mixture (KOH & KzCO3) is dissolved in distilled water completely. It is titrated against 0.11 HClcaqs. An average of 120 cm3 HCl is used. Find KOH % by mass.

KOH + HC| 
$$\rightarrow$$
 KC| + H<sub>2</sub>O  $\Rightarrow$  KOH : HC| = 1:1  
K<sub>2</sub>CO<sub>3</sub> + 2HC|  $\rightarrow$  2KC|<sub>2</sub> + CO<sub>2</sub> + H<sub>2</sub>O  $\Rightarrow$  K<sub>2</sub>CO<sub>3</sub> : HC| = 1:2  
Let x g be KOH mmss, (0.75- $\chi$ ) g be K<sub>2</sub>CO<sub>3</sub> mmss.  
HC| in KOH va + HC| in K<sub>2</sub>CO<sub>3</sub> va = 0.1 × 0.12  

$$\frac{\chi}{39.1+16+1} + \frac{0.75 - \chi}{39.1\cdot2+12+16\cdot3} \cdot 2 = 0.012$$

$$\frac{\chi}{56.1} + \frac{0.75 - \chi}{69.1} = 0.012$$

$$\frac{69.1\chi}{3876.51} + \frac{0.75 - \chi}{69.1} = 0.012$$

$$13\chi = 46.51812-42.075$$

$$\chi = 0.342$$

KOH% by mass =  $\frac{0.342}{0.75} \cdot 100\% = 45.6\%$